


Weighing Bench for laboratory tables.

A heavy, highly durable Granite Stone Plate that is polished and smooth.
The Shock-absorbing Tape under the Plate minimises vibration and also prevents the undermining of chemicals.

The polished and sealed surface also allows to use the Weighing Bench in clean rooms.
The leveling is done by the feet of the established scales or gadgetries.




Design:

Art. L58745001 Weighing Bench for laboratory tables.

Bevelled edge, thickness 3 cm, surface polished
weight : ca. 25 KG, dimensions: 54,0 cm x 54,0 cm

Art. L58745002 Weighing Bench for laboratory tables.

Bevelled edge, thickness 3 cm, surface polished
weight : ca. 26 KG, dimensions: 65,0 cm x 50,0 cm



Equipment, Optional:

Art. L58745005 Rack for Weighing Bench

This flat stainless steel rack comes with damping pedestals.

The Weighing Bench (54,0 cm x 54,0 cm) has to be inserted in the rack.

By placing the assembly on your laboratory table you will get a superb weighing table.