

EBI 11 Starter Kit SK 6100

FOR TEMPERATURE AND PROCESS MONITORING

NEW

EBI 11 Data Logger Starter Kit SK 6100

1,586.00 €*

*valid until end of 2026



SK 6100 consists of:

- 1x EBI 11-T210
- 1x Software Winlog.pro
- 1x Interface IF 300

EBI 11 Data Logger Starter Kit SK 6100

The **Starter Kit SK 6100** contains all necessary parts of the ebro portfolio to get started with temperature and process monitoring. The **high-quality and versatile** design of the **EBI 11-T210**, as well as the high precision and fast response time make the logger ideal for the evaluation of sterilization and pasteurization processes. The **Winlog.pro evaluation software** enables the integration of several data loggers as well as the evaluation and visualization of the recorded data. Thanks to the four readout ports of the **IF 300 interface**, a subsequent increase in the number of loggers is easily possible.

Measurement range	-30 °C ... +140 °C (-22 °F ... + 284 °F)
Accuracy	± 0.1 °C (0 °C ... +140 °C)
Sampling rate	1 s ... 24 hours
Data memory	15,000 measurement values
3-point factory calibration certificate	0 °C, +60 °C and +134 °C
External temperature probe	20 mm length, Ø 3 mm, axial pointed, M5 Thread, Pt 1000

Type	Description	Part No.	Discount
SK 6100	1 x EBI 11-T210 data logger + 1x Winlog.pro + 1x IF 300	1250-6100	C01

Headquarter:

Xylem Analytics Germany Sales GmbH & Co. KG, ebro
Peringerstr. 10 · 85055 Ingolstadt Germany
Phone +49 841 95478-0 · Fax +49 841 95478-80
Internet: www.ebro.com · E-Mail: ebro@xylem.com

Representatives worldwide:

www.ebro.com

ebro® is a registered trademark of Xylem Inc. or one of its subsidiaries. Technical changes reserved
© 2022 Xylem Analytics Germany Sales GmbH & Co. KG.

December 2025

-ebro-
a xylem brand